

# AbsoluteIDQ<sup>®</sup> Stero17 kit requirements

This document specifies additional lab equipment and chemicals that are required but <u>not</u> included with the kit. Before starting with the kit, ensure that all required items are available.

## LC-MS system

SCIEX	Waters			
Mass spectrometer				
SCIEX 5500 or 6500+ series (QTRAP or Triple Quad) with TurboV™ ion source	Waters Xevo TQ-XS, TQ-S or TQ-S micro with ESI source			
LC and autosampler				
<ul> <li>Standard flow HPLC or UHPLC system with column oven.</li> <li>NanoLC and microLC are not supported!</li> <li>96-well plate autosampler with temperature control (10 °C)</li> <li>Injection volume range: 10 – 20 µL</li> <li>Special sample rack for deep-well plates required for the following autosamplers:</li> <li>Shimadzu: part no. 228-37546-92 (Deep-well MTP Rack, DWP-96)</li> <li>UltiMate 3000: Thermo part no. 6820.4083 (30 – 36 mm height)</li> </ul>	<ul> <li>Waters ACQUITY UPLC<sup>®</sup> system with sample manager, solvent manager, and column oven:</li> <li>Classic (loop-based)</li> <li>I-Class (FTN)</li> <li>H-Class (FTN)</li> <li>Injection volume range: 5 – 10 µL</li> </ul>			
Tubing				
<ul> <li>Recommended between injector and column system: steel capillary</li> <li>Recommended between column system and mass spectrometer: red PEEK Tubing 1/16" OD x 0.005" ID</li> </ul>				

## LC-MS instrument status

#### Condition

- Serviced, calibrated and tuned instruments according to manufacturer's recommendations on a regular basis
- LC-MS system must be free of ion-pair reagent

## Solvents and chemicals

Solvents and chemicals	Purity	
Acetonitrile, isopropanol, methanol, water	LC-MS grade	
Dichloromethane	p.a. grade of higher	



## Analytical column

HPLC		
Analytical column	AbsoluteIDQ <sup>®</sup> Stero17 kit HPLC column: biocrates product number 9120052120301	
Precolumn	AbsoluteIDQ <sup>®</sup> Stero17 kit HPLC precolumn: biocrates product number 9120052120349	
Precolumn holder	AbsoluteIDQ <sup>®</sup> Stero17 kit HPLC setup box: biocrates product number 9120052120318	

UHPLC		
Analytical column	AbsoluteIDQ <sup>®</sup> Stero17 kit UHPLC column: biocrates product number 9120052120356	
Precolumn	AbsoluteIDQ <sup>®</sup> Stero17 kit UHPLC precolumn: biocrates product number 9120052120363	
Precolumn holder	recolumn holder AbsoluteIDQ® Stero17 kit UHPLC setup box: biocrates product number 9120052120370	

## Software

SCIEX	Waters	
MS software		
Analyst <sup>®</sup> version 1.5.2 or later	MassLynx® version 4.1 or later	
Data processing computer		

#### Data processing computer

- Oracle® database:
  - Option 1 (recommended): full commercial version, updates and support provided by Oracle. Expert required.
  - Option 2: Oracle Database Express Edition (Oracle XE), free of charge, provided by biocrates. No updates and support provided by Oracle.
- MetIDQ<sup>™</sup> software (provided by biocrates)
- Windows 7 or later, 64-bit architecture
- Microsoft .NET Framework 4.5.2 or later
- Microsoft Visual C++ redistributables (32 & 64 bit each): 2005, 2008, 2010, 2012, and 2013



# Laboratory equipment

	Laboratory equipment			
Shaker	Any model with adjustable speed (450 – 1200 rpm) including a tray for plates and vials. Recommended: Eppendorf MixMate® or ThermoMixer®			
Pipettes	<ul> <li>Repeater, e.g. Eppendorf Multipette<sup>®</sup> E3 (with 1.0 and 10 mL tips) or similar electronic model</li> <li>Single channel: volume range 10 μL – 1000 μL</li> <li>8-Channel: 1000 μL</li> </ul>			
Vortexer	Any model			
Ultrasonic bath	Any model			
Solvent bottles	50 – 1000 mL			
Drying and extraction – Option 1				
Nitrogen evaporator for 96-well plates*		Vacuum-based system		
<ul> <li>Examples:</li> <li>Porvair blowdown evaporator MiniVap®</li> <li>Techne (FSC496D)</li> <li>Sample concentrator from BenchTop Lab Systems (96 samples, BT1604)</li> <li>VLM evaporators</li> <li>Organomation MICROVAP microplate evaporator</li> <li>Biotage TurboVap® 96</li> </ul>		<ul> <li>Vacuum pump, e.g.: <ul> <li>Waters WAT085114, WAT085123 or WAT085115</li> </ul> </li> <li>Vacuum manifold for 96-well SPE plates, e.g.: <ul> <li>Extraction plate manifold for Oasis 96 Well Plates, Waters 186001831</li> <li>PlatePrep 96-well vacuum manifold, starter kit, Sigma-Aldrich 575650-U</li> </ul> </li> <li>nd installation in fume hood</li> </ul>		
	Requires nitrogen supply and installation in fume hood Drying and extraction – Option 2			
Pressure manifold	for 96-well plates*	Adapter for pressure manifold		
<ul> <li>Examples:</li> <li>Waters positive pressure-96 processor (186006961)</li> <li>Biotage<sup>®</sup> PRESSURE+ 96 manifold (PPM-96)</li> <li>CEREX<sup>®</sup> system 96 processor (288-0001)</li> <li>Agilent positive pressure manifold 96 processor (PPM-96)</li> <li>TECAN Resolvex<sup>®</sup> M10 96</li> </ul>		biocrates product number 9120052120981		
Requires nitrogen supply and installation in fume hood				

\* See document "Technical guide-Nitrogen evaporators and pressure manifolds" for details.