

## MxP<sup>®</sup> Quant 1000 kit requirements

This document specifies additional lab equipment and chemicals that are required but **not** included with the kit. Before starting with the kit, ensure that all required items are available.

### LC-MS system

Agilent	SCIEX	Waters
<b>Mass spectrometer</b>		
<b>6495C</b> Triple Quadrupole LC/MS System	<b>5500+, 6500+ or 7500</b> series (QTRAP or Triple Quad) with TurboV™ ion source	<b>Xevo TQ-XS</b> with ESI source
<b>LC system and autosampler</b>		
<ul style="list-style-type: none"> <li>– Agilent 1290 Infinity I/II UHPLC system with column oven.</li> <li>– 800 bar injection valve (Agilent part no. 5067-6698) replaces standard 1300 bar valve</li> <li>– 1290 Infinity II inline filters, 0.3 µm (2 pcs., before and after injection valve, Agilent part no. 5067-6189)</li> <li>– 96-well plate autosampler with temperature control (4 °C)</li> <li>– Injection volume range: 2 – 20 µL</li> </ul>	<ul style="list-style-type: none"> <li>– Standard flow UHPLC system with column oven.</li> <li>– 96-well plate autosampler with temperature control (4 °C)</li> <li>– Injection volume range: 1–20 µL</li> <li>– Shimadzu and ExionLC autosamplers capable of accepting 96 well plates: <a href="#">Plate calibration video</a></li> </ul>	<ul style="list-style-type: none"> <li>– Waters ACQUITY UPLC<sup>®</sup> system with sample manager, solvent manager, and column oven:               <ul style="list-style-type: none"> <li>– Classic (loop-based)</li> <li>– I-Class (FTN)</li> <li>– H-Class (FTN)</li> </ul> </li> <li>– Injection volume range: 1–20 µL</li> <li><b>FTN samplers require extension loop of 50 µL or needle size of 30 µL!</b></li> </ul>
<b>NanoLC and microLC are not supported!</b>		
<b>Tubing</b>		
<ul style="list-style-type: none"> <li>– Recommended between injector and analytical column: steel capillary</li> <li>– Recommended between analytical column and mass spectrometer: red PEEK Tubing 1/16" OD x 0.005" ID</li> </ul>		
<b>Condition</b>		
<ul style="list-style-type: none"> <li>– <b>Instrument is well maintained and serviced</b> on a regular basis according to manufacturer's recommendations</li> <li>– <b>MS tuned and calibrated in both positive <u>and</u> negative ion modes</b></li> <li>– LC-MS system must be <b>free of ion pairing reagents</b></li> </ul>		

## Column components (available as package or individually)

Component	Requirement
Column package	MxP® Quant 1000 kit column (1 pc.) + precolumn (3 pcs.) biocrates part number 22558
Analytical column	MxP® Quant 1000 kit column biocrates part number 22565
Precolumn	MxP® Quant 1000 kit precolumn biocrates part number 22572

## Laboratory equipment

Equipment	Requirements (may be available from local distributor)	
	Nitrogen evaporator for 96-well plates	Pressure manifold for 96-well plates
Nitrogen evaporator* <b>or</b> pressure manifold* – Requires nitrogen supply – Must be in a fume hood – No heater required	<b>Examples:</b> – Porvair blowdown evaporator MiniVap® – Techne (FSC496D) – Sample concentrator from BenchTop Lab Systems (96 samples, BT1604) – VLM evaporators – Organomation MICROVAP microplate evaporator – Biotage TurboVap® 96	<b>Examples:</b> – Porvair UltraPPM LITE – Waters positive pressure-96 processor (186006961) – Biotage® PRESSURE+ 96 manifold (PPM-96) – CEREX® system 96 processor (288-0001) – Agilent positive pressure manifold 96 processor (PPM-96) – TECAN Resolvex® M10 96 <div style="background-color: #fff9c4; padding: 5px; border: 1px solid #ccc;">  <b>Adapter required for drying plate biocrates part number 20981</b> </div>
Centrifuge	Must be able to centrifuge 96-well plates of 5 cm height at 500 x g	Not required when a pressure manifold is used.
Plate heater	For 96-well plate incubation at 35 °C, e.g. Eppendorf shaker with heated lid (ThermoTop®) or oven	
Shaker	Any model with adjustable speed (450–1200 rpm) including a tray for plates and vials. Recommended: Eppendorf MixMate® or ThermoMixer®	
Pipettes	– Repeater, e.g. Eppendorf Multipipette® E3 (with 2.5 and 10 mL tips) or similar electronic model – Single channel: volume range 10 µL–1000 µL – 8-Channel: 10 µL–150 µL	
Vortexer	Any model	
Balance	Accuracy < 1 mg	
Solvent bottles	50–1000 mL	

\* See document “Technical guide-Nitrogen evaporators and pressure manifolds” for details.

## Solvents and chemicals

Solvents and chemicals	Purity
Methanol, water, acetonitrile, isopropanol	LC-MS grade
Formic acid	LC-MS grade (e.g. Honeywell Fluka™ 56302-50 mL) <b>Fresh! New bottle or opened within past 6 months</b>
Phenyl isothiocyanate (PITC)	99%, for protein sequencing (e.g. Sigma-Aldrich 317861) <b>Fresh! New bottle or opened within past 6 months</b>
Pyridine	>99% (e.g. Acros Organics 131780500, 50 mL) <b>Fresh! New bottle or opened within past 6 months</b>
Triethylamine	>99% (e.g. Sigma-Aldrich 90335-100 mL) <b>Fresh! New, unopened bottle</b>
Ammonium acetate	LC-MS grade
Phosphate buffered saline (PBS)	p.a. grade (e.g. Sigma-Aldrich P4417)

## Software

Agilent	SCIEX	Waters
MS software		
<b>MassHunter</b> version 10.1 or later	<b>Analyst®</b> version 1.7.3 or later <b>SCIEX OS</b> version 3.4 or later	<b>MassLynx®</b> version 4.2 or later
biocrates WebIDQ workflow manager		

Subscriptions available through biocrates webshop or local distributor:

- |  |                               |
|--|-------------------------------|
| – <b>WebIDQ basic</b> subscription (4 GB storage space)  | biocrates part number 22336.1 |
| – <b>WebIDQ core</b> subscription (8 GB storage space)   | biocrates part number 21629.3 |
| – <b>WebIDQ core+</b> subscription (50 GB storage space) | biocrates part number 21926.3 |
| – <b>WebIDQ on-premises</b> subscription                 | biocrates part number 21636.3 |

For detailed system requirements please visit:

<https://shop.biocrates.com/faq/What-are-the-system-requirements-for-WebIDQ>