



biocrates  
The future of research and health

# SMartIDQ alpha kit

Efficient workflow for busy laboratories

[biocrates.com](https://biocrates.com)

## Mass spectrometry-based metabolomics kits



High efficiency



Quantitative, reproducible



Streamlined integration



Comprehensive analysis

For research use only | not for use in diagnostic procedures

# SMartIDQ alpha kit

The SMartIDQ alpha kit helps you streamline your lab operations with meaningful metabolite coverage and high-throughput capabilities



## High efficiency and reproducibility

- ▶ Minimal sample volumes (only 10 µL)
- ▶ ONE ready-to-use kit for multiple matrices
- ▶ Little personnel requirement
- ▶ Automatable sample preparation



## Comprehensive analysis

- ▶ Set of well-characterized metabolites
- ▶ Discover biological insights with MetaboINDICATOR™ – 138 metabolite sums and ratios
- ▶ Leading quantitative metabolite database



## Coverage


- ▶ 107 metabolites from 13 biochemical classes

## Sample volume | type

- ▶ 10 µl | For use with various sample matrices and species

## Throughput

- ▶ 480 samples/week/instrument



## Kit components

- ▶ Patented 96-well filter plates with integrated internal standards
- ▶ Calibration standards and quality controls
- ▶ Test samples for system suitability test
- ▶ Quick-start guide and user manual
- ▶ Protocol for different matrices
- ▶ Instrument-specific methods

[shop.biocrates.com/kit-products/](https://shop.biocrates.com/kit-products/)

## Read more



biocrates' standardized  
metabolomics technology

Get more information on how to generate  
reproducible metabolomics data

[biocrates.com/our-technology](https://biocrates.com/our-technology)



SMartIDQ alpha kit

Streamlined lab operations with meaningful  
metabolite coverage

[biocrates.com/smartidq-alpha-kit](https://biocrates.com/smartidq-alpha-kit)