

AbsoluteIDQ[®] Stero17 kit

Ready-to-use, standardized solution for the quantitation of 17 steroid hormones.

[biocrates.com](https://www.biocrates.com)

Mass spectrometry-based metabolomics kits



Ready-to-use



Automated workflow



Standardized, quantitative



High throughput

AbsoluteIDQ® Stero17 kit



Comprehensive coverage in steroid analysis

- ▶ Quantitatively measures 17 steroid hormones
- ▶ Covers all five steroid hormone classes
- ▶ The benchmark technology for targeted steroid hormone screening



Superior data quality

- ▶ Excellent accuracy and precision
- ▶ Proven in proficiency tests

Kit components



Reagents & consumables

- ▶ 96-well SPE plate
- ▶ Internal & calibration standards
- ▶ Quality controls
- ▶ System suitability sample

Methods & protocols

- ▶ Instrument-specific acquisition and quantification methods
- ▶ Manual with detailed descriptions

Workflow manager MetIDQ™

- ▶ Project design
- ▶ Automated technical validation and data normalization
- ▶ Basic statistics (optional)

Read more



biocrates' standardized metabolomics technology

Get more information on how to produce reproducible metabolomics data

biocrates.com/our-technology

Last revised 03/ 2021



Applications of metabolomics

Learn how metabolomics is already changing lives

biocrates.com/applications